

# Pressure Mapping

Pressure mapping is a specialised measurement technology used to measure and visualize the contact pressure distribution between the human body and a supporting surface and equipment interface, e.g. person, chair or sling. Care & Independence commission independent pressure mapping experts to conduct such trials to ascertain sling performance and help identify areas of risk. The subsequent scientific data insight has enabled Care & Independence to develop solutions and vastly improve upon the areas which indicate tissue viability risks, pain or other health concerns to the equipment user.



## RESULTS FOR GLOVE™ AIRFLOW (GAIR)

As with all our pressure mapping test results, the Polyair® cushion is the better performing cushion when a GLOVE™ in-chair style is added into the mix. However the high air permeability and pressure-reducing fabric benefits of the GLOVE™ Airflow sling is clearly evident in the test results for the Integrity® cushion too.

In the baseline mapping tests where the subject was sat clothed direct upon both cushion types, a clear low pressure result was returned as indicated by the expanse of blue colouring. [Fig. 1 & 2]

In the two minute seating test [Fig. 3 & 4], the pressure performance on both cushions with the addition of the GLOVE™ Airflow sling return results that are either **the same or better average** than the baseline test results.

A ten minute seating test returns similar results [Fig. 5 & 6] but with the Polyair®/GAIR combination showing a **further enhanced pressure reduction** on both two minute and baseline results.



MEDICAL CUSHION TYPE:

1. Polyair® comfort cushion
2. Sumed Integrity® Static High Risk

SUBJECT: Male, 5'6", 82.5kg

DATE OF TEST: April 2021

TESTER: Sumed International (UK) Ltd

\*mmHg stands for millimetres of mercury and is used as a pressure measurement.